

Anti-ITGAM antibody (1-100) [S5134RM] (STJ11105134)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Rabbit
Reactivity	Mouse/Rat

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:1000 ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 127kDa Observed Mw: 170kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3684
Gene Symbol	ITGAM
UniProt ID	ITAM_HUMAN
Immunogen Region	1-100
Immunogen Sequence	MALRVLTLTALTLCHEFNLD TENAMTFQENARGFGQSVVQ LQGSRVVVGAPQEIVAANQR GSLYQCDYSTGSCPEIRLQV PVEAVNMSLGLSLAATTSP
Specificity	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human CD11b/ITGAM (NP_000623.2).

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Western blot analysis of lysates from Mouse spleen, using CD11b/ITGAM Rabbit monoclonal antibody (STJ11105134) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 μ g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 180s.