

Anti-PIM1 antibody (213-313) (STJ11105343)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:1000 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 35kDa/45kDa Observed Mw: 44kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5292
Gene Symbol	PIM1
UniProt ID	PIM1_HUMAN
Immunogen Region	213-313
Immunogen Sequence	IRYHRYHGRSAAVWSLIGILL YDMVCGDIPFEHDEEIIIRGQ VFFRQVRSSECQHLIRWCLA LRPSDRPTFEEIQNHPWMQD VLLPQETAEIHLHSLSPGPS K
Specificity	A synthetic peptide corresponding to a sequence within amino acids 213-313 of human PIM1 (NP_002639.1).

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot analysis of Rat kidney, using PIM1 Rabbit polyclonal antibody (STJ11105343) at 1:600 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25µg per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 60s.