

## Anti-CD317 antibody (55-152) (STJ11105374)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Rabbit
Reactivity	Mouse

### PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:1000-1:5000 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 19kDa Observed Mw: 28-40kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">69550</a>
Gene Symbol	<a href="#">Bst2</a>
UniProt ID	<a href="#">BST2_MOUSE</a>
Immunogen Region	55-152
Immunogen Sequence	SVACRDGLRAQAECRNTHL LQRQLTRTQDSLLQAETQAN SCNLTVVTLQESLEKKVSQA LEQQARIKELENEVTKLNQE LENLRIQKETSSTVQVNS
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 55-152 of mouse CD317 (NP_932763.1).

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot analysis of Mouse ovary, using CD317 Rabbit polyclonal antibody (STJ11105374) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25µg per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 30s.