

Anti-DCAF1 antibody (1058-1396) (STJ11105384)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:1000 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 118kDa/168kDa/169kDa Observed Mw: 170kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	9730
Gene Symbol	DCAF1
UniProt ID	DCAF1_HUMAN
Immunogen Region	1058-1396
Immunogen Sequence	ASFPKYGGVDGGCFDRHLIF SRFRPISVVFREANEDESGFT CCAFSARERFLMLGTCTGQL KLYNVFSGQEEASYNCHNSA ITHLEPSRDGSLLLTSATWS QPLSALWGMKSVFDMKHSFT EDHYVEFSKHSQDRVIGTKG DIAHIYDIQTGNKLLTLFNP DLANNYKRNCATFNPTDDL LNDGVLWDVRSQAIIHKFDK FNMNISGVFHPNGLEVIINT EIWDLRTFHLLHTVPALDQ
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1058-1396 of human VPRBP (NP_055518.1).

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Western blot analysis of various lysates, using VPRBP Rabbit polyclonal antibody (STJ11105384) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 180s.