



Anti-FOXQ1 antibody (1-100) (STJ11105413)

GENERAL INFORMATION	
Product Type	Primary antibodies
Applications	IHC-P/ELISA
Host / Source	Rabbit
Reactivity	Human
PRODUCT PROPERTIES	
Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IHC-P:1:50-1:200 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 41kDa Observed Mw: Refer to figures
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
TARGET INFORMATION	
Gene ID	94234
Gene Symbol	FOXQ1
UniProt ID	FOXQ1_HUMAN
Immunogen Region	1-100
Immunogen Sequence	MKLEVFVPRAAHGDKQGSDDL EGAGGSDAPSPLSAAGDDSL GSDGDCAANSPAAGGGGARDT QGDGEQSAGGGPGAEEAIPA AAAAAVVAEGAEAGAAGPGA
Specificity	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human FOXQ1 (NP_150285.3).
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Immunohistochemistry analysis of paraffin-embedded human lung cancer using FOXQ1 Rabbit polyclonal antibody (STJ11105413) at dilution of 1:50 (40x lens). Perform high pressure antigen retrieval with 10 mM citrate buffer pH 6.0 before commencing with immunohistochemistry staining protocol.



Immunohistochemistry analysis of paraffin-embedded human cervical using FOXQ1 Rabbit polyclonal antibody (STJ11105413) at dilution of 1:50 (40x lens). Perform high pressure antigen retrieval with 10 mM citrate buffer pH 6.0 before commencing with immunohistochemistry staining protocol.