

Anti-PRELID3B antibody (80-180) (STJ11105451)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:1000 ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 21kDa Observed Mw: 21kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	51012
Gene Symbol	PRELID3B
UniProt ID	PLD3B_HUMAN
Immunogen Region	80-180
Immunogen Sequence	EHSVVDPEKTMELKSTNIS FTNMVSVDERLIYKPHPQDP EKT VLTQEAITVKG VSLSS YLEGLMASTISSNASKGREA MEWVIHKLNAEIEELTASAR G
Specificity	A synthetic peptide corresponding to a sequence within amino acids 80-180 of human PRELID3B (NP_057129.2).

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot analysis of lysates from 293T cells using PRELID3B Rabbit polyclonal antibody (STJ11105451) at 1:1000 dilution. Secondary antibody:HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 μ g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection:ECL Basic Kit. Exposure time:60s.