

Anti-ACSM3 antibody (28-130) (STJ11105568)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:1000 ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 66kDa Observed Mw: 66kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6296
Gene Symbol	ACSM3
UniProt ID	ACSM3_HUMAN
Immunogen Region	28-130
Immunogen Sequence	LHKDNRTATPQNFSNYESMK QDFKLGIPYFNFAKDVLQDQ WTDKEKAGKKPSNPAFWWIN RNGEEMRWSFEELGSLSRKF ANILSEACSLQRGDRVILIL PRV
Specificity	A synthetic peptide corresponding to a sequence within amino acids 28-130 of human ACSM3 (NP_005613.2).

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot analysis of lysates from Mouse kidney using ACSM3 Rabbit polyclonal antibody (STJ11105568) at 1:1000 dilution. Secondary antibody:HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 μ g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection:ECL Basic Kit. Exposure time:90s.