

Anti-Cy3/Cy5 antibody [S5637RM] (STJ11105637)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	DB/ELISA
Host / Source	Rabbit
Reactivity	Species Independent

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	DB:1:1000-1:5000 ELISA: Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Immunogen Sequence	Cy5-NHS-KLH
Specificity	Chemical compounds corresponding to Cy5.

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Dot-blot analysis of all sorts of compounds using Cy3/Cy5 Rabbit monoclonal antibody (STJ11105637) at 1:2000 dilution. NA NA Lane 1 : Cy5 Goat Anti-Mouse IgG (H+L) NA NA Lane 2 : Cy5-NHS-BSA NA NA Lane 3 : Cy3 Goat Anti-Mouse IgG (H+L) NA NA Lane 4 : Cy3-NHS-BSA NA NA Lane 5 : FITC Goat Anti-Mouse IgG (H+L) NA NA Lane 6 : HRP Goat Anti-Mouse IgG (H+L)