

Anti-NECTIN3 antibody (127-366) (STJ11105695)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:100-1:500 ELISA: Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 61kDa Observed Mw: 83kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	25945
Gene Symbol	NECTIN3
UniProt ID	NECT3_HUMAN
Immunogen Region	127-366
Immunogen Sequence	SLNDATITLHNIGFSDSGKY ICKAVTFPLGNAQSSTTVTV LVEPTVSLIKGPDSLIDGNG ETVAAICIAATGKPVAHIDW EGDLGEMESTTTSFPNETAT IISQYKLFPTRFARGRRITC VVKHPALEKDIRYSFILDQ YAPEVSVTG YDGNWVFVGRKG VNLKCNADANPPPFKSVWSR LDGQWPDGLLASDNTLHFVH PLTFNYSGVYICKVTNSLGQ RSDQKVYISAYNSVASLN
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 127-366 of human PVRL3 (NP_001230215.1).

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Western blot analysis of Mouse testis, using PVRL3 Rabbit polyclonal antibody (STJ11105695) at 1:400 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25µg per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 45s.