

Anti-PTPRB antibody (1173-1260) (STJ11105723)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:1000 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 224kDa Observed Mw: 224kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5787
Gene Symbol	PTPRB
UniProt ID	PTPRB_HUMAN
Immunogen Region	1173-1260
Immunogen Sequence	PASVSHLRGSNRNTTDSLWF NWSPASGDFDFYELILYNPN GTKKENWKDKDLTEWRFQGL VPGRKYVLWVVTHSGDLSNK VTAESRTA
Specificity	Recombinant Protein corresponding to a sequence within amino acids 1173-1260 of human PTPRB (NP_002828.3).

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot analysis of lysates from U-251 MG cells using PTPRB Rabbit polyclonal antibody (STJ11105723) at 1:800 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 Mu g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.