

Anti-ANGPTL1 antibody (1-100 aa) (STJ11106376)
GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 57kDa Observed Mw: Refer to figures
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	9068
Gene Symbol	ANGPTL1
UniProt ID	ANGL1_HUMAN
Immunogen Region	1-100 aa
Immunogen Sequence	MKTFWTWLGVLFFLLVDTGH CRGGQFKIKKINQRRYPRAT DGKEEAKKCAYTFLVPEQRI TGPICVNTKGQDASTIKDMI TRMDLENLKDVLRSQKREID
Specificity	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human ANGPTL1 (NP_004664.1).

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---