

Anti-ANKRD1 antibody (1-130 aa) [SRM5089] (STJ11106872)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:1000 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 36kDa Observed Mw: 43kDa

TARGET INFORMATION

Gene ID	27063
Gene Symbol	ANKRD1
UniProt ID	ANKR1_HUMAN
Immunogen Region	1-130 aa
Immunogen Sequence	MMVLKVEELVTGKKNGNGEA GEFLPEDFRDGEYEA AVTLE KQEDLKTLLAHPVTLGEQQW KSEKQREAE LKKKKLEQRSK LENLEDLEIIQLKKRKKYR KTKVPVVKEPEPEIITEPVD VPTFLKAALE
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-130 of human ANKRD1 Rabbit PolymAb® (NP_055206.2).

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---