

Anti-HADHA antibody (545-763 aa) [SRM5146] (STJ11106984)

GENERAL INFORMATION	
Product Type	Primary antibodies
Applications	WB/IHC-P/IF/ICC/IP/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat
PRODUCT PROPERTIES	
Clonality	Monoclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:1000 IHC-P:1:50-1:200 IF/ICC:1:50-1:200 IP:0.5 Mu g-4 Mu g antibody for 200 Mu g-400 Mu g extracts of whole cells ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specif
Formulation	PBS with 0.05% Proclin300, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 83kDa Observed Mw:
TARGET INFORMATION	
Gene ID	3030
Gene Symbol	HADHA
UniProt ID	ECHA_HUMAN
Immunogen Region	545-763 aa
Immunogen Sequence	FYTTRCLAPMMSEVIRILQE GVDPKKLDLSLTTSFGFPVGA ATLVDEVGVDVAKHVAEDLG KVFGERFGGGNPELLTQMVS KGFLGRKSGKGFYIYQEGVK RKDLNSDMDASILASLKLPPK SEVSSDEDIQFRLVTRFVNE AVMCLQEGILATPAEGDIGA VFGLGFPPCLGGPFRFVDLY GAQKIVDRLKKYEAAYGKQF TPCQLLADHANSPNKKFYQ
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 545-763 of human HADHA (NP_000173.2).
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.