

Anti-L1R antibody (1-152 aa) (STJ11107073)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Rabbit
Reactivity	Monkeypox

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:1000 ELISA: Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 18kDa Observed Mw: 17kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Immunogen Region	1-152 aa
Immunogen Sequence	MDHNQYLLTMFFADDDSFYK YFASQDDESSLSDILOITQY LDFLLLLLIQSKNKLEAVGH CYESLSEYRQLTKFTDSQD FKFLFNKVPITDGRVKLNK GYLDFVISLMRFKKESALA TTAIDPVRYIDPRRDIAFSN VMDILKSNKVEQ
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-152 of Monkeypox virus Zaire-96-I-16 L1R (NP_536512.1).

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Western blot analysis of Recombinant Monkeypox virus L1R Protein (STJP001528) using L1R Rabbit pAb (STJ11107073) at 1:1000 dilution.
 Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution.
 Lysates/proteins: 10ng/5ng/1ng per lane.
 Blocking buffer: 3% nonfat dry milk in TBST.
 Detection: ECL Basic Kit
 Exposure time: 1s.