

Anti-Pan-Lactic Acid-Lysine antibody [SRM5229] {Magnetic Beads} (STJ11107111)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	IP
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat/Other

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	Lot specific
Conjugation	Magnetic Beads
Purification	Affinity purification
Dilution Range	IP:30ul antibody (bead slurry) for 200 Mu g-400 Mu g extracts of whole cells
Formulation	PBS with 0.05% Proclin300, pH 7.3.
Isotype	IgG
Molecular Weight	Refer to figures

TARGET INFORMATION

Specificity	A synthetic peptide corresponding to a sequence containing Pan Lactic acid-Lysine.
-------------	--

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---