

**Anti-ADSL antibody (1-484 aa) [SRM5256] (STJ11107158)**
**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Rabbit
Reactivity	Human

**PRODUCT PROPERTIES**

Clonality	Monoclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:1000 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 55kDa Observed Mw: Refer to figures
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

**TARGET INFORMATION**

Gene ID	<a href="#">158</a>
Gene Symbol	<a href="#">ADSL</a>
UniProt ID	<a href="#">PUR8_HUMAN</a>
Immunogen Region	1-484 aa
Immunogen Sequence	MAAGGDHGGSPDSYRSPLASR YASPEMCFVFSDDRYKFRTWR QLWLWLAEEQTLGLPITDE QIQEMKSNLENIDFKMAAEE EKRLRHDMAHVHTFGHCCP KAAGIIHLGATSCYVGDNTD LIILRNALDLLLPKLARVIS RLADFAKERASLPTLGFTHF QPAQLTTVGKRCCLWIQDLC MDLQNLKRVRRDRLFRGVKG TTGTQASFLQLFEGDDHKVE QLDKVMTEKAGFKRAFIIT
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-484 of human ADSL (NP_000017.1).

**ADDITIONAL INFORMATION**

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---