

Anti-MARCH5 antibody (179-278 aa) [SRM5286] (STJ11107188)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	WB/IHC-P/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:1000 IHC-P:1:50-1:200 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 31kDa Observed Mw: Refer to figures
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	54708
Gene Symbol	MARCH5
UniProt ID	MARH5_HUMAN
Immunogen Region	179-278 aa
Immunogen Sequence	LNSIFPGIGCPVPRIPAEAN PLADHVSATRILCGALVFPT IATIVGKLMFSSVNSNLQRT ILGGIAFVAIKGAFKVYFKQ QQYLRQAHRKILNYPEQEEA
Specificity	A synthetic peptide corresponding to a sequence within amino acids 179-278 of human MARCH5 (NP_060294.1).

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---