

Anti-MYH9 antibody (1861-1960 aa) [SRM5291] (STJ11107193)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:1000 ELISA: Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 227kDa Observed Mw: Refer to figures
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4627
Gene Symbol	MYH9
UniProt ID	MYH9_HUMAN
Immunogen Region	1861-1960 aa
Immunogen Sequence	YKDQADKASTRLKQLKRQLE EAEEEEAQRANASRRKLQREL EDATETADAMNREVSSLKNK LRRGDLFPVVPRRMARKGAG DGSDEEVDGKADGAEAKPAE
Specificity	A synthetic peptide corresponding to a sequence within amino acids 1861-1960 of human MYH9 (NP_002464.1).

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---