

Anti-ICAM-1/CD54 antibody [SRM5305] (STJ11107208)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	WB/IHC-P/IF/ICC/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:1000-1:30000 IHC-P:1:2000-1:10000 IF/ICC:1:2000-1:4000 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.09% Sodium Azide, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 58kDa/59kDa Observed Mw: 82-100kDa

TARGET INFORMATION

Gene ID	15894 3383
Gene Symbol	Icam1 ICAM1
UniProt ID	ICAM1_MOUSE ICAM1_HUMAN
Specificity	Recombinant protein of human/mouse ICAM-1/CD54.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---