

Anti-CD66a/CEACAM1 antibody [SRM5306] (STJ11107209)

GENERAL INFORMATION	
Product Type	Primary antibodies
Applications	WB/IHC-P/IF/ICC/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse
PRODUCT PROPERTIES	
Clonality	Monoclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:3000-1:12000 IHC-P:1:4000-1:10000 IF/ICC:1:2000-1:8000 ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.09% Sodium Azide, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 30kDa/37kDa/43kDa/46kDa/50kDa/57kDa Observed Mw: 100-180kDa
TARGET INFORMATION	
Gene ID	26365 634
Gene Symbol	Ceacam1 CEACAM1
UniProt ID	CEAM1_MOUSE CEAM1_HUMAN
Specificity	Recombinant protein of human/mouse CD66a/CEACAM1.
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.