

Anti-TAZ antibody (1-270 aa) [SRM5317] (STJ11107221)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC-P/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:1500-1:9000 IHC-P:1:50-1:200 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with Proclin300, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 44kDa Observed Mw: 60kDa

TARGET INFORMATION

Gene ID	25937
Gene Symbol	WWTR1
UniProt ID	WWTR1_HUMAN
Immunogen Region	1-270 aa
Immunogen Sequence	MNPASAPPLPPPGQQVIHV TQDLTDLEALFNSVMNPKP SSWRKKILPESFFKEPDSGS HSRQSSTDSSGGHPGPRLAG GAQHVRSHSSPASLQLGTGA GAAGSPAQQHAHLRQQSYDV TDELPLPPGWEMTFTATGQR YFLNHIEKITTWQDPRKAMN QPLNHMNLHPAVSSTPVPQR SMAVSQPNLVMNHQHQQQMA PSTLSQQNHPTQNPPAGLMS MPNALTTQQQQQKLRQLR
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-270 of human TAZ (NP_056287.1)

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---