

Anti-PSMA/FOLH1 antibody (50-150 aa) [SRM5319] (STJ11107223)

GENERAL INFORMATION	
Product Type	Primary antibodies
Applications	WB/IHC-P/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat
PRODUCT PROPERTIES	
Clonality	Monoclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:3000-1:6000 IHC-P:1:100-1:500 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.09% Sodium Azide, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 84kDa Observed Mw: 100kDa
TARGET INFORMATION	
Gene ID	2346
Gene Symbol	FOLH1
UniProt ID	FOLH1_HUMAN
Immunogen Region	50-150 aa
Immunogen Sequence	TNITPKHNMKAFLDELKAEN IKKFLYNFTQIPHLAGTEQN FQLAKQIQSQWKEFGLDSVE LAHYDVLLSYPNKTHPNYIS IINEDGNEIFNTSLFEPPPP G
Specificity	A synthetic peptide corresponding to a sequence within amino acids 50-150 of human PSMA/FOLH1 (NP_004467.1).
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.