

Anti-FGF1 antibody (56-155 aa) [SRM5328] (STJ11107234)

GENERAL INFORMATION	
Product Type	Primary antibodies
Applications	WB/IHC-P/IF/ICC/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat
PRODUCT PROPERTIES	
Clonality	Monoclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:1000-1:30000 IHC-P:1:2000-1:20000 IF/ICC:1:200-1:800 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.09% Sodium Azide, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 17kDa Observed Mw: 17kDa
TARGET INFORMATION	
Gene ID	2246
Gene Symbol	FGF1
UniProt ID	FGF1_HUMAN
Immunogen Region	56-155 aa
Immunogen Sequence	HIQLQLSAESVGEVYIKSTE TGQYLAMDTDGLLYGSQTPN EECLFLERLEENHYNTYISK KHAENWVFLKKNKGSKRG PRTHYGQKAILFLPLPVSSD
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 56-155 of human FGF1 (P05230).
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.