

Anti-ADAM9 antibody (29-210 aa) [SRM5357] (STJ11107266)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:6500-1:13000 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.09% Sodium Azide, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 91kDa Observed Mw: 100kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	8754
Gene Symbol	ADAM9
UniProt ID	ADAM9_HUMAN
Immunogen Region	29-210 aa
Immunogen Sequence	AARPGFQQTSHLSSYEIITP WRLTRERREAPRPYSKQVSY VIQAEGKEHIIHLERNKDLL PEDFVVYTYNKEGTLITDHP NIQNHCHYRGYVEGVHNSSI ALSDCFGLRGLLHLENASYG IEPLQNSSHFEHIIYRMDDV YKEPLKCGVSNKDIEKETAK DEEEPPSMTQLLRRRRRAVL PQ
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 29-210 of human ADAM9 (NP_003807.1).

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---