

## Anti-PFKP antibody (336-540 aa) [SRM5361] (STJ11107271)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

### PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:1000-1:4000 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.09% Sodium Azide, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 86kDa Observed Mw: 86kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">5214</a>
Gene Symbol	<a href="#">PFKP</a>
UniProt ID	<a href="#">PFKAP_HUMAN</a>
Immunogen Region	336-540 aa
Immunogen Sequence	ATPDTPACVVSLNGNHAVRL PLMECVQMTQDVQKAMDERR FQDAVRLRGRSFAGNLNTYK RLAIKLPDDQIPKTNVAV INVGAPAAGMNAAVRSVAVRV GIADGHRMLAIYDGFDFGFAK GQIKEIGWTDVGGWTGQGGG ILGTKRVLPGKYLEEATQM RTHSINALLIIGFEAYLGL LELSAAREKHEEFCVPMVMV PATVS
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 336-540 of human PFKP (NP_002618.1).

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---