

Anti-IP6K1 antibody (1-220 aa) [SRM5369] (STJ11107279)
GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:2500-1:10000 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.09% Sodium Azide, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 50kDa Observed Mw: 60kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	9807
Gene Symbol	IP6K1
UniProt ID	IP6K1_HUMAN
Immunogen Region	1-220 aa
Immunogen Sequence	MCVCQTMEVQGQYGNASRAG DRGVLLPEFPIHQVGGHSSMM RYDDHTVCKPLISREQRFYE SLPPEMKEFTPEYKGVVSVCFEGDSDGYINLVAYPYVESE TVEQDDTTEREQPRRKHSRR SLHRSGSGSDHKEEKASLSLETSESSQEAQSPKVELHSHSEVPPFQMLDGNSSGLSSEKISH NPWSLRCHKQQLSRMRSESK DRKLYKFLLENNVHHFKYP
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-220 of human IP6K1 (NP_695005.1).

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---