

Anti-PSMC4 antibody (1-165 aa) [SRM5381] (STJ11107292)
GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-P/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IHC-P:1:200-1:800 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.09% Sodium Azide, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 47kDa Observed Mw: Refer to figures

TARGET INFORMATION

Gene ID	5704
Gene Symbol	PSMC4
UniProt ID	PRS6B_HUMAN
Immunogen Region	1-165 aa
Immunogen Sequence	MEEIGILVEKAQDEIPALSV SRPQTGLSFLGPEPEDLEDL YSRYKKLQQELEFLEVQEEY IKDEQKNLKKFLHAQEEVK RIQSIPLVIGQFLEAVDQNT AIVGSTTGSNYVYRILSTID RELLKPNASVALHKHSNALV DVLPEADSSIMMLTSDQKP DVMYA
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-165 of human PSMC4 (NP_006494.1).

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---