

Anti-RXRB antibody (1-100 aa) [SRM5409] (STJ11107320)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	IHC-P/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IHC-P:1:200-1:800 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.09% Sodium Azide, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 57kDa Observed Mw: Refer to figures

TARGET INFORMATION

Gene ID	6257
Gene Symbol	RXRB
UniProt ID	RXRB_HUMAN
Immunogen Region	1-100 aa
Immunogen Sequence	MSWAARPPFLPQRHAAGQCG PVGVRKEMHCGVASRWRRRR PWLDPAAAAAAVAGGEQQT PEPEPGEAGRDGMGDSGRDS RSPDSSSPNPLPQGVPPSP
Specificity	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human RXRB (NP_068811.1).

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---