

Anti-GMNN antibody (1-209 aa) [SRM5439] (STJ11107352)

GENERAL INFORMATION	
Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat
PRODUCT PROPERTIES	
Clonality	Monoclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:1500-1:3000 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.09% Sodium Azide, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 24kDa Observed Mw: 30kDa
TARGET INFORMATION	
Gene ID	51053
Gene Symbol	GMNN
UniProt ID	GEMI_HUMAN
Immunogen Region	1-209 aa
Immunogen Sequence	MNPSMKQKQEEIKENIKNSS VPRRTLKMIQPSASGSLVGR ENELSAGLSKRKHRNDHLTS TTSSPGVIVPESSKNLGG VTQESFDLMIKENPSSQYWK EVAEKRRKALYEALKENEKL HKEIEQKDNEIARLKKENKE LAEVAEHVQYMAELIERLNG EPLDNFESLDNQEFDSEET VEDSLVEDSEIGTCAEGTVS SSTDAKPCI
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-209 of human GMNN (NP_056979.1).
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.