

## Anti-LRP6 antibody (200-300 aa) (STJ11107425)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Rabbit
Reactivity	Mouse/Rat

### PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:1000 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 180kDa Observed Mw: 160kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">4040</a>
Gene Symbol	<a href="#">LRP6</a>
UniProt ID	<a href="#">LRP6_HUMAN</a>
Immunogen Region	200-300 aa
Immunogen Sequence	DAKLNFIHKSNLDGTRQAV VKGSLPHPFALTLFEDILYW TDWSTHSILACNKYTGEGLR EIHSDIFSPMDIHAFSQQRQ PNATNPCGIDNGGCSHLCLM S
Specificity	A synthetic peptide corresponding to a sequence within amino acids 200-300 of human LRP6 (NP_002327.2).

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot analysis of various lysates using LRP6 Rabbit pAb (STJ11107425) at 1:1000 dilution.

Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution.

Lysates/proteins: 25 Mu g per lane.

Blocking buffer: 3% nonfat dry milk in TBST.

Detection: ECL Basic Kit

Exposure time: 90s.