

Anti-Mouse IgG antibody [SMM5490] (STJ11107555)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IP
Host / Source	Mouse
Reactivity	Mouse

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Protein A/G purification
Dilution Range	IP:0.5ug-4ug antibody for 200ug-400ug extracts of whole cells
Formulation	PBS with 0.09% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Specificity	None
-------------	------

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Immunoprecipitation of over-expressed Myc-tagged protein in 293T cells incubated using Myc-tag antibody (STJ11103700). A mock served as negative control using mouse Control IgG (STJ11107555) and over-expressed 293T cell lysate served as positive control.