

Anti-Rabbit IgG antibody {Magnetic Beads} (STJ11107561)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IP
Host / Source	Rabbit
Reactivity	Rabbit

PRODUCT PROPERTIES


Concentration	Lot specific
Conjugation	Magnetic Beads
Purification	Protein A/G purification
Dilution Range	IP:30ul antibody (bead slurry) for 200 Mu g-400 Mu g extracts of whole cells
Formulation	PBS with 0.05% Proclin300, pH 7.3.
Isotype	IgG
Storage Instruction	Store at 4°C for up to 1 year from the date of receipt. Do not freeze.

TARGET INFORMATION

Specificity	None
-------------	------

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Immunoprecipitation analysis of 300 ug extract cell lysate from 293T cells transfected with GSK3B expression vector containing a C-terminal DDDDK tag (1x) with 30uL Magnetic Beads-conjugated Rabbit anti-DDDDK-tag mAb antibody (STJ11107584). Magnetic Beads-conjugated Rabbit IgG isotype control mAb was used as a negative control (STJ11107561). Western blot was performed from the immunoprecipitate using Mouse anti DDDDK-Tag mAb antibody (STJ11103697) at 1:2000 dilution. Secondary antibody: HRP-conjugated AffiniPure Goat Anti-Mouse IgG Light Chain (STJS001198).