

## Anti-MBP-Tag antibody [SMM5511] {HRP} (STJ11107581)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB
Host / Source	Mouse
Reactivity	Species Independent

### PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	Lot specific
Conjugation	HRP
Purification	Affinity purification
Dilution Range	WB:1:2000-1:10000
Formulation	PBS with 50% Glycerol, pH 7.3.
Isotype	IgG1k
Molecular Weight	Protein Mw: 42.5 kDa Observed Mw: Refer to Figures
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Specificity	A synthetic peptide of MBP-Tag
-------------	--------------------------------

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---