

Anti-MYC-Tag antibody [SRM5518] (STJ11107591)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IP
Host / Source	Rabbit
Reactivity	Species Independent

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IP:1:50-1:200

TARGET INFORMATION

Specificity	A synthetic peptide corresponding to Myc-Tag.
-------------	---

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---