

## Anti-Integrin alpha 5 (ITGA5/CD49e) antibody [SRM13055] (STJ11107785)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC-P/FC/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

### PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB, 1:1000-1:6000 IHC-P, 1:200-1:2000 FC (intra) , 1:50-1:200 ELISA, Recommended starting concentration is 1 $\mu$ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 0.05% BSA, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">3678</a>
Gene Symbol	<a href="#">ITGA5</a>
UniProt ID	<a href="#">ITA5_HUMAN</a>
Immunogen Sequence	Human

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	-------------------------------------------------------------------------------------------------------------------------