

Anti-MonoMethyl-Histone H3-K79 antibody [SRM13169] (STJ11107899)

GENERAL INFORMATION	
Product Type	Primary antibodies
Applications	WB/DB/IHC-P/ELISA/ChIP/ChIP-seq
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat/Other
PRODUCT PROPERTIES	
Clonality	Monoclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB, 1:4000-1:16000 DB, 1:2000-1:20000 IHC-P, 1:1000-1:5000 ELISA, Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements. ChIP, 5 µg antibody for 5 µg-10 µg o
Formulation	PBS containing 50% Glycerol and 0.05% BSA, preserved with proclin300 or Sodium Azide (as specified on the Certificate of Analysis) , pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
TARGET INFORMATION	
Gene ID	8290 8350/8351/8352/8353/8354/8355/8356/8357/8358/8968
Gene Symbol	H3-4 H3C1.H3C2.H3C3.H3C4.H3C6.
UniProt ID	H31T_HUMAN H31_HUMAN
Immunogen Sequence	Human
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.