

Anti-SMURF1 antibody [SMM13262] (STJ11107992)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB, 1:500-1:1000 ELISA, Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.1% Sodium Azide and 0.1% gelatin.
Isotype	IgG2ak
Storage Instruction	Store at 4°C for up to 1 year from the date of receipt. Do not freeze.

TARGET INFORMATION

Gene ID	57154
Gene Symbol	SMURF1
UniProt ID	SMUF1_HUMAN
Immunogen Sequence	Human

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---