

## Anti-HEXA antibody [SRM13268] (STJ11107998)

### GENERAL INFORMATION

|               |                     |
|---------------|---------------------|
| Product Type  | Primary antibodies  |
| Applications  | WB/IF-P/IHC-P/ELISA |
| Host / Source | Rabbit              |
| Reactivity    | Human/Mouse/Rat     |

### PRODUCT PROPERTIES

|                     |                                                                                                                                                                                                    |
|---------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Clonality           | Monoclonal                                                                                                                                                                                         |
| Concentration       | Lot specific                                                                                                                                                                                       |
| Conjugation         | Unconjugated                                                                                                                                                                                       |
| Purification        | Affinity purification                                                                                                                                                                              |
| Dilution Range      | WB, 1:500-1:1000<br>IF-P, 1:50-1:200<br>IHC-P, 1:1000-1:5000<br>ELISA, Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements. |
| Formulation         | PBS containing 50% Glycerol and 0.05% BSA, preserved with proclin300 or Sodium Azide (as specified on the Certificate of Analysis) , pH 7.3.                                                       |
| Isotype             | IgG                                                                                                                                                                                                |
| Storage Instruction | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.                                                                                                     |

### TARGET INFORMATION

|                    |                            |
|--------------------|----------------------------|
| Gene ID            | <a href="#">3073</a>       |
| Gene Symbol        | <a href="#">HEXA</a>       |
| UniProt ID         | <a href="#">HEXA_HUMAN</a> |
| Immunogen Sequence | Human                      |

### ADDITIONAL INFORMATION

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|