

**Anti-CD352/SLAMF6 antibody [SRM13343] (STJ11108073)****GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	WB/IHC-P/FC/ELISA
Host / Source	Rabbit
Reactivity	Human

**PRODUCT PROPERTIES**

Clonality	Monoclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB, 1:1000-1:4000 IHC-P, 1:200-1:800 FC, 1:500-1:1000 ELISA, Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS containing 50% Glycerol and 0.05% BSA, preserved with proclin300 or Sodium Azide (as specified on the Certificate of Analysis) , pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

**TARGET INFORMATION**

Gene ID	<a href="#">114836</a>
Gene Symbol	<a href="#">SLAMF6</a>
UniProt ID	<a href="#">SLAF6_HUMAN</a>
Immunogen Sequence	Human

**ADDITIONAL INFORMATION**

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---