

## Anti-CD3E antibody [SRM13366] (STJ11108096)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC-P/mlIHC/FC
Host / Source	Rabbit
Reactivity	Human

### PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB, 1:1000-1:5000 IHC-P, 1:50-1:200 mlIHC, 1:50-1:200 FC, 1:100-1:500
Formulation	PBS with 0.09% Sodium Azide, 0.05% BSA, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">916</a>
Gene Symbol	<a href="#">CD3E</a>
UniProt ID	<a href="#">CD3E_HUMAN</a>
Immunogen Sequence	Human

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	-------------------------------------------------------------------------------------------------------------------------