

## Anti-IL12B antibody [SRM13385] (STJ11108115)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | WB/IHC-P/FC/ELISA  |
| Host / Source | Rabbit             |
| Reactivity    | Human/Mouse        |

### PRODUCT PROPERTIES

|                     |                                                                                                                                                                                                  |
|---------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Clonality           | Monoclonal                                                                                                                                                                                       |
| Concentration       | Lot specific                                                                                                                                                                                     |
| Conjugation         | Unconjugated                                                                                                                                                                                     |
| Purification        | Affinity purification                                                                                                                                                                            |
| Dilution Range      | WB, 1:2000-1:8000<br>IHC-P, 1:200-1:800<br>FC, 1:100-1:500<br>ELISA, Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements. |
| Formulation         | PBS containing 50% Glycerol and 0.05% BSA, preserved with proclin300 or Sodium Azide (as specified on the Certificate of Analysis) , pH 7.3.                                                     |
| Isotype             | IgG                                                                                                                                                                                              |
| Storage Instruction | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.                                                                                                   |

### TARGET INFORMATION

|                    |                             |
|--------------------|-----------------------------|
| Gene ID            | <a href="#">16160</a>       |
| Gene Symbol        | <a href="#">IL12b</a>       |
| UniProt ID         | <a href="#">IL12B_MOUSE</a> |
| Immunogen Sequence | Mouse                       |

### ADDITIONAL INFORMATION

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|