

Anti-CD4 antibody [SRM13418] (STJ11108148)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IP/FC/ELISA
Host / Source	Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB, 1:500-1:1000 IP, 0.5 μ g-4 μ g antibody for 400 μ g-600 μ g extracts of whole cells FC, 1:100-1:500 ELISA, Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirement
Formulation	PBS with 0.03% proclin300, 0.2% BSA, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	12504
Gene Symbol	Cd4
UniProt ID	CD4_MOUSE
Immunogen Sequence	Mouse

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---