

## Anti-CD45RB antibody [SRM13440] (STJ11108170)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IF/ICC/IF-P/IHC-P/FC/ELISA
Host / Source	Rabbit
Reactivity	Human

### PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB, 1:2000-1:8000 IF/CC, 1:50-1:200 IF-P, 1:50-1:200 IHC-P, 1:50-1:200 FC, 1:500-1:1000 ELISA, Recommended starting concentration is 1 $\mu$ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS containing 50% Glycerol and 0.05% BSA, preserved with proclin300 or Sodium Azide (as specified on the Certificate of Analysis) , pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">5788</a>
Gene Symbol	<a href="#">PTPRC</a>
UniProt ID	<a href="#">PTPRC_HUMAN</a>
Immunogen Sequence	Human

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---