

**Anti-CD162/PSGL-1 antibody [SRM13471] (STJ11108201)****GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	ELISA/WB/IF/ICC/FC
Host / Source	Rabbit
Reactivity	Human

**PRODUCT PROPERTIES**

Clonality	Monoclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB, 1:30000-1:120000 IF/ICC, 1:100-1:1600 FC, 1:50-1:200
Formulation	PBS containing 50% Glycerol and 0.05% BSA, preserved with proclin300 or Sodium Azide (as specified on the Certificate of Analysis) , pH 7.3.
Isotype	IgG
Storage Instruction	Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

**TARGET INFORMATION**

Gene ID	<a href="#">6404</a>
Gene Symbol	<a href="#">SEPLG</a>
UniProt ID	<a href="#">SELPL_HUMAN</a>
Immunogen Sequence	Human

**ADDITIONAL INFORMATION**

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---