

Anti-ICAM-1/CD54 antibody [SRM13497] (STJ11108227)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC-P/FC/ELISA
Host / Source	Rabbit
Reactivity	Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB, 1:5000-1:20000 IHC-P, 1:200-1:800 FC, 1:500-1:1000 ELISA, Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS containing 50% Glycerol and 0.05% BSA, preserved with proclin300 or Sodium Azide (as specified on the Certificate of Analysis) , pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	25464
Gene Symbol	Icam1
UniProt ID	ICAM1_RAT
Immunogen Sequence	Rat

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---