

Anti-CD162/PSGL-1 antibody [SRM13502] (STJ11108232)
GENERAL INFORMATION

| | |
|---------------|---------------------|
| Product Type | Primary antibodies |
| Applications | IF-P/IHC-P/FC/ELISA |
| Host / Source | Rabbit |
| Reactivity | Mouse/Rat |

PRODUCT PROPERTIES

| | |
|---------------------|---|
| Clonality | Monoclonal |
| Concentration | Lot specific |
| Conjugation | Unconjugated |
| Purification | Affinity purification |
| Dilution Range | IF-P, 1:200-1:2000 IHC-P, 1:200-1:2000 FC, 1:100-1:500 ELISA, Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements. |
| Formulation | PBS containing 50% Glycerol and 0.05% BSA, preserved with proclin300 or Sodium Azide (as specified on the Certificate of Analysis) , pH 7.3. |
| Isotype | IgG |
| Storage Instruction | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles. |

TARGET INFORMATION

| | |
|--------------------|-------|
| Immunogen Sequence | Mouse |
|--------------------|-------|

ADDITIONAL INFORMATION

| | |
|------|---|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|---|