

Anti-Phospho-Uroplakin III/UPK3A-S10 antibody (STJ11108322)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IF/ICC/IHC-P/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB, 1:10000-1:50000 IF/CC, 1:300-1:3000 IHC-P, 1:3000-1:12000 ELISA, Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.09% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	22270
Gene Symbol	Upk3a
UniProt ID	UPK3A_MOUSE
Immunogen Sequence	Mouse

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---