

Anti-Phospho-Histone H3-S10 antibody [SRM13680] (STJ11108410)

GENERAL INFORMATION		
	Product Type	Primary antibodies
	Applications	WB/IHC-P/IF/ICC/ELISA
	Host / Source	Rabbit
	Reactivity	Human/Mouse/Rat/Other
PRODUCT PROPERTIES		
	Clonality	Monoclonal
	Concentration	Lot specific
	Conjugation	Unconjugated
	Purification	Affinity purification
	Dilution Range	WB, 1:500-1:1000 IHC-P, 1:50-1:200 IF/ICC, 1:50-1:200 ELISA, Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
	Formulation	PBS with 0.02% Sodium Azide, 0.05% BSA, 50% Glycerol, pH7.3.
	Isotype	IgG
	Storage Instruction	Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
TARGET INFORMATION		
	Gene ID	8290 8350/8351/8352/8353/8354/8355/8356/8357/8358/8968
	Gene Symbol	H3-4 H3C1.H3C2.H3C3.H3C4.H3C6.
	UniProt ID	H31T_HUMAN H31_HUMAN
	Immunogen Sequence	Human
ADDITIONAL INFORMATION		
	Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.