

**Anti-Phospho-GSK3 Beta + GSK3 Alpha-Y216+Y279 antibody [SRM13708] (STJ11108438)**
**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Rabbit
Reactivity	Mouse

**PRODUCT PROPERTIES**

Clonality	Monoclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB, 1:1000-1:4000 ELISA, Recommended starting concentration is 1 $\mu$ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.09% Sodium Azide, 0.05% BSA, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

**TARGET INFORMATION**

Gene ID	<a href="#">2932</a> <a href="#">2931</a>
Gene Symbol	<a href="#">GSK3B</a> <a href="#">GSK3A</a>
UniProt ID	<a href="#">GSK3B_HUMAN</a> <a href="#">GSK3A_HUMAN</a>
Immunogen Sequence	Human

**ADDITIONAL INFORMATION**

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---