

## Anti-CD40 antibody {ABflo® 488} (STJ11108444)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Rabbit
Reactivity	Human

### PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	ABflo® 488
Purification	Affinity purification
Dilution Range	FC, 5 Mu l per 10 <sup>6</sup> cells in 100 Mu l volume
Formulation	PBS containing 0.2% BSA, preserved with proclin300 or sodium azide (as specified on the Certificate of Analysis) , pH 7.3.
Isotype	IgG
Molecular Weight	31kDa
Storage Instruction	Store at 4°C for up to 1 year from the date of receipt. Do not freeze.

### TARGET INFORMATION

Gene ID	<a href="#">958</a>
Gene Symbol	<a href="#">CD40</a>
UniProt ID	<a href="#">TNRS_HUMAN</a>

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	-------------------------------------------------------------------------------------------------------------------------